


DECON-CYCLE[®]

**Sterile Pharmaceutical Clean Room Formula
Low pH Phenolic Germicidal Detergent**



SEE THE NEW
Veltek
SimpleMix[®]
SYSTEM



Veltek Associates, Inc.

15 Lee Boulevard • Malvern, PA 19355-1234

Tel: (610)644-8335 Fax: (610)644-8336 E-mail: vai@sterile.com VL-501-98 12/18/03

DECON-CYCLE®

Sterile Pharmaceutical Clean Room Formula

PRODUCT DESCRIPTION

DECON-CYCLE LOW pH Phenolic is a phosphate free, dilutable, hospital detergent effective in 10 minutes at 20 degrees Celsius, in hard water up to 400 ppm (calculated as CaCO₃), in the presence of 5% blood serum, used to simulate organic soil. DECON-CYCLE is recommended for use by pharmaceutical and biotechnology industries medical device manufacturers hospitals and any health care institutions that are dedicated to controlling the hazards of cross contamination in their establishments. DECON-CYCLE's multi-phenolic formula is designed to clean, disinfect and deodorize any washable inanimate surface in one easy step. When used as directed, DECON-CYCLE is strong enough to kill a broad spectrum of pathogenic bacteria, yet mild enough to have no harmful effect on the surface being disinfected.

DECON-CYCLE is produced sterile by filtration at 0.2 microns in a GMP Class 100 manufacturing area and subsequently gamma radiated to a 10⁻⁶ SAL level. DECON-CYCLE is available sterile and non-sterile in 2 ounce bottles and a 1 gallon containers. Each container is double bag packaged and delivered each time with the lot specific Certificate of Analysis and Sterility Report. DECON-CYCLE is completely traceable and has been completely validated for sterility and shelf life.

DECON-CYCLE is a phosphate free germicidal detergent that has efficacy performance which has been confirmed using A.O.A.C. protocol testing at a use dilution of 1:256 in the presence of 5% blood serum (used to simulate organic soil load) at 20 degrees Celsius with an exposure time of 10 minutes.

The SimpleMix® System Technology Alternative:

As a innovating step to address this problematic situation, Veltek Associates, Inc has developed the SimpleMix® System Technology. The system provide for the transfer of the sterile disinfectant and sterile water in a sealed container to the aseptic area. The system container is double bag packaged for easy transfer and eliminated all internal and external concerns for the use of a sterile, compliant disinfectant or sporicide.

The Patent Pending SimpleMix® System provides a sealed multi-chamber container that when activated mixes the two solutions. The top part contains the concentrate disinfectant or sporicide and the bottom part contains the USP WFI Quality Water. **Just open the cap, pull the tab and they instantly mix together.** The system design permits the easy transfer of the sterile disinfectant or sporicide and sterile water to the aseptic manufacturing area without concern for the transfer of contamination See **VAI Technical Document (SMix-01 Rev. 02/03)**

Uses: For use on hard, inanimate surfaces in aseptic filling and gowning rooms, general manufacturing areas, or on machinery, tables, counters, laminar flow benches, floors, walls, stainless steel, porcelain, glass and chrome.

ORDERING INFORMATION

<u>Order#</u>	<u>Description</u>	<u>Quan/cs.</u>
DCY-01	DECON-CYCLE 1 Gallon Non-Sterile	4
DCY-02	DECON-CYCLE 1 Gallon Sterile	4
DCY-03-2Z	DECON-CYCLE 2 oz. Sterile	24
DCY-03-4Z	DECON-CYCLE 4 oz. Sterile	24
DCY-04-1/2Z	DECON-PHENE SimpleMix Sterile 1 gallon	4
DCY-05-1/2Z	DECON-PHENE SimpleMix Non-Sterile 1 gallon	4

Available Technical Data Supplements (Upon Request)

SC-001: Sterile Chemical Product Guide

VL-501: DECON-CYCLE Technical Data Report

Veltek Associates, Inc.

15 Lee Boulevard • Malvern, PA 19355-1234

Tel: (610)644-8335 Fax: (610)644-8336 E-mail: vai@sterile.com VL-501-98 12/18/03

DECON-CYCLE[®]

Sterile Pharmaceutical Clean Room Formula
Low pH Phenolic Germicidal Detergent

PRODUCT LABELING

Veltek Associates, Inc.

15 Lee Boulevard • Malvern, PA 19355-1234

Tel: (610)644-8335 Fax: (610)644-8336 E-mail: vai@sterile.com VL-501-98 12/18/03

O V E R V I E W

4

(CENTER PANEL)

DECON-CYCLE®

Sterile Pharmaceutical Clean Room Formula

Low pH Phenolic

Germicidal Detergent Effective In Hard Water
up to 400 Ppm (Calculated as CaCO₃)

With Organic Soil Tolerance

Tuberculocidal - Germicidal - Pseudomonacidal
Fungicidal - Virucidal

Formulated for Effective Aseptic Controlled Environment (Clean Room) Sanitation

ACTIVE INGREDIENTS:

Ortho-Phenylphenol.....	8.085%
Orlhobenzyl-Benzyl-para-Chlorophenol.....	6.650%
INERT INGREDIENTS.....	85.265%

SEE LEFT PANEL FOR ADDITIONAL PRECAUTIONARY STATEMENTS

Keep Out of Reach of Children

DANGER

Veltek Associates, Inc.
15 Lee Boulevard
Malvern, PA 19355-1234
Tel: (610)644-8335

EPA Reg. No. 211-62-068959
EPA Est. No. 068959-PA-01

Veltek Associates, Inc.
15 Lee Boulevard • Malvern, PA 19355-1234
Tel: (610)644-8335 Fax: (610)644-8336 E-mail: vai@sterile.com VL-501-98 12/18/03

O V E R V I E W

5

(LEFT PANEL)

PRECAUTIONARY STATEMENTS **Hazards to Humans and Domestic Animals**

DANGER

Keep out of reach of children. Corrosive. Causes eye and skin damage. Do not get in eyes, on skin, or on clothing. Avoid breathing spray mist. Harmful if swallowed. Wear goggles or face shield and rubber gloves when handling. Avoid contamination of food. Remove and wash contaminated clothing,

STATEMENT OF PRACTICAL TREATMENT

In case of eye contact, Immediately flush eyes or skin with plenty of water for at least 15 minutes. For eyes, call a physician. Remove and wash all contaminated clothing before reuse. If swallowed, drink large quantities of water. Call a physician.

NOTE TO PHYSICIAN

Probable mucosal damage may contraindicate the use of gastric lavage. Measures against shock, respiratory depression, and convulsion may be needed,

ENVIRONMENTAL HAZARDS

Do not discharge effluent containing this product into lakes, streams, ponds, estuaries, oceans, or other waters unless in accordance with the requirements of a National Pollutant Discharge Elimination System (NPDES) permit and the permitting authority has been notified in writing prior to discharge. Do not discharge effluent containing this product to sewer systems without previously notifying the local sewage treatment plant authority. For guidance contact your State Water Board or Regional Office of the EPA. This product is toxic to fish and aquatic organisms.

STORAGE AND DISPOSAL

Do not contaminate water, food, or food by storage or disposal.

STORAGE - Store in a securely closed original container. Store in area that will prevent cross contamination of other pesticides, fertilizer, food, and food. Avoid storage at temperature extremes or in sunlight.

PESTICIDE DISPOSAL - Pesticide wastes are acutely hazardous. Improper disposal of excess pesticide, spray mixture, or rinsate is a violation of Federal Law. If these wastes cannot be disposed of by use according to label instructions, contact your state Pesticide or Environmental Control Agency, or the hazardous waste representative at the nearest EPA Regional office for guidance.

CONTAINER DISPOSAL - Triple rinse (or equivalent), then offer for recycling or reconditioning, or puncture and dispose of in a sanitary landfill, by incineration, or, if allowed by state and local authorities, by burning. If burned, stay out of smoke.

Veltek Associates, Inc.

15 Lee Boulevard • Malvern, PA 19355-1234

Tel: (610)644-8335 Fax: (610)644-8336 E-mail: vai@sterile.com VL-501-98 12/18/03

O V E R V I E W

6

(REQUIRED)

DIRECTION FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with it's labeling

This product is not to be used as a terminal sterilant / high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

CLEANING AND DISINFECTION

For broad spectrum disinfection of gram negative (as demonstrated by efficacy against Salmonella cholerasuis) and gram positive (as demonstrated by efficacy against Staphylococcus aureus and Pseudomonas aeruginosa) bacteria, add 1/2 ounce DECON-CYCLE per 1 gallon of water. Remove heavy soil or gross filth then apply solution with a mop, cloth, sponge, or hand-pump trigger sprayer so as to wet surface thoroughly. Allow to remain wet for 10 minutes and then air dry. If higher detergency is desired, increase dilution to 1 ounce per gallon of water.

To clean and disinfect toilet bowl and urinal surfaces, remove gross soils from surfaces with bowl brush and flush toilet or urinal prior to applying the use solution. Remove residual bowl water from toilet bowl and urinal by forcing water over the trap with a bowl swab. Press bowl swab against side of the bowl or urinal to remove excess water from the applicator. Add 1/2 ounce of DECON-CYCLE directly into the bowl or urinal trap water. Swab or brush the interior surface of the fixture thoroughly, making sure to get under the rim. For disinfection, allow to stand for 10 minutes then flush. Wipe exterior surfaces with cloth or sponge.

(HEALTH CARE FACILITY USE)

DECON-CYCLE is a phosphate free, dilutable, hospital germicidal-detergent effective in 10 minutes at 20° C, in hard water up to 400 ppm (calculated as CaCO₃), in the presence of 5% blood serum, used to simulate organic soil. DECON-CYCLE's multi-phenolic formula is designed to clean, disinfect, and deodorize any washable, inanimate, non-porous surface in one easy step. When used as directed, DECON-CYCLE is strong enough to kill a broad spectrum of pathogenic bacteria, yet It is mild enough to have no harmful effects on the surface being disinfected, when used as directed.

DECON-CYCLE has been formulated for use on hard, nonporous, inanimate surfaces commonly found in facilities and areas such as: hospitals, nursing homes, other health care institutions, and establishments that are dedicated to controlling the hazards of cross contamination. This product is specifically designed for use in hospitals, veterinary clinics, animal shelters, kennels, hospital emergency departments, ambulance and patient transfer vehicles, police and law enforcement vehicles, prisoner holding areas, offices, funeral homes, restrooms, medical and dental offices, and day care facilities.

DECON-CYCLE is recommended for use on hard, inanimate surfaces such as door handles, medical bed surfaces, springs, wheelchairs, walls, floors, light switches, linen carts, stretcher wheels, hampers, telephones, clean up carts, bassinets, dressing carts, toilet bowl surfaces, urinal surfaces, showers, bathtubs, bed frames, mattress springs, impervious vinyl surfaces, counter tops, table tops, sinks, waste containers, lockers, desks, patient transfer lift, patient scales, examination tables, dental chairs, medical lamps, stethoscopes, bedpans, blood donor chairs, crutches, instrument trays, operating tables, and other hard, non-porous surfaces that need disinfection.

When used as directed, DECON-CYCLE may be used on hard, nonporous surfaces, materials, and components such as:

Veltek Associates, Inc.

15 Lee Boulevard • Malvern, PA 19355-1234

Tel: (610)644-8335 Fax: (610)644-8336 E-mail: vai@sterile.com VL-501-98 12/18/03

O V E R V I E W

7

Metals: Stainless Steel, Chrome Plated Steel, Brass, copper and Aluminum
Plastics: Fiberglass™, Formica™, Nylon, ABS Lexan™, Plexiglass™, Acrylics and Vinyl
Porcelain, Glass and Ceramic Tile
Latex, Enamel, & Epoxy Painted Surfaces
Vinyl, Rubber, Terrazzo, Marble, Quarry Tile Floor Surfaces
Urethane, Lacquer, and other surfaces not harmed by water based cleaners.

DECON-CYCLE is a dilutable germicidal detergent formulated for use where broad spectrum disinfection is required in the presence of moderate amounts of organic soil, in hard water up to 400 ppm (calculated as CaCO₃).

EFFICACY PERFORMANCE

The broad spectrum efficacy of DECON-CYCLE has been confirmed using A.O.A.C. protocol testing at a use-dilution of 1:256 in hard water up to 400 ppm (calculated as CaCO₃) in the presence of 5% blood serum (used to simulate organic soil load), at 20° C with an exposure time of 10 minutes.

GERMICIDAL:

DECON-CYCLE is effective against a broad spectrum of Gram positive and Gram negative pathogens as confirmed by the A.O.A.C. Use-Dilution test, in hard water up to 400 ppm (calculated as CaCO₃), and in the presence of 5% organic soil load at 20° C with an exposure time of 10 minutes against the following organisms:

Staphylococcus aureus	Escherichia coli
Pseudomonas aeruginosa	Salmonella choleraesuis
Methicillin Resistant Staphylococcus aureus (MRSA)	

FUNGICIDAL:

DECON-CYCLE is effective against Trichophyton mentagrophytes pathogenic fungi (the athletes foot fungus), as confirmed by the A.O.A.C. Use-Dilution test, in hard water up to 400 ppm (calculated as CaCO₃), and in the presence of 5% organic soil load at 20° C with an exposure time of 10 minutes, in use areas such as shower room floors, locker room benches, bathmats, or any other hard, nonporous, inanimate environmental surfaces.

TUBERCULOCIDAL:

DECON-CYCLE is effective against Mycobacterium bovis (BCG) as demonstrated through testing using the A.O.A.C. Method "Confirmative in Vitro Test for Determining Tuberculocidal Activity", in hard water up to 400 ppm (calculated as CaCO₃) and in the presence of 5% organic soil load, at 20° C with an exposure time of 10 minutes.

VIRUCIDAL*

DECON-CYCLE is effective against the following viruses on hard, nonporous, inanimate surfaces, as confirmed by the A.O.A.C. use-dilution test, in hard water up to 400 ppm (calculated as CaCO₃) and in the presence of 5% organic soil load, at 20° C with an exposure time of 10 minutes:

HIV-1 (AIDS virus)
Herpes Simplex, Type 2

Veltek Associates, Inc.

15 Lee Boulevard • Malvern, PA 19355-1234

Tel: (610)644-8335 Fax: (610)644-8336 E-mail: vai@sterile.com VL-501-98 12/18/03

O V E R V I E W
8

***VIRUCIDAL**

***KILLS HIV-1 ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type-I (HIV-1) (associated with AIDS).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 (AIDS VIRUS) OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS.

PERSONAL PROTECTION: Clean-up should always be done wearing protective latex gloves, gowns, masks, and eye protection.

CLEANING PROCEDURE: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of the disinfectant.

CONTACT TIME: Leave surface wet for 10 minutes with a 1/2 ounce per gallon use-dilution (1:256), DECON-CYCLE and water

DISPOSAL OF INFECTIOUS MATERIALS: Blood, body fluids, cleaning materials, and clothing should be autoclaved and disposed according to local regulations for infectious waste disposal.

PELIGRO: Si no puede leer en ingles, pregunte a su supervisor sobre las instrucciones de uso apropiadas antes de trabajar con este producto.

(Alternative label text for center panel)

**Statements to be added only to label of containers
destined for aseptic controlled environments (clean rooms)**

Product and container sterilized and distributed by _____

Sterile Pharmaceutical Clean Room Formula

DECON-CYCLE is acceptable for use by pharmaceutical and biotechnology industries medical device manufacturers hospitals and any health care institutions that are dedicated to controlling the hazards of cross contamination in their establishments. For use on hard, inanimate surfaces in aseptic filling and gowning rooms, general manufacturing areas, or on machinery, tables, counters, laminar flow benches, floors, walls, stainless steel, porcelain, glass, and chrome.

FUNGICIDAL: At 1/2 ounce per gallon of water, DECON-CYCLE is also effective against Trichophyton mentagrophytes pathogenic fungi in use areas as aseptic filling and gowning rooms, and general manufacturing areas or on any other hard, non-porous, inanimate environmental surfaces.

This product is not to be used as a terminal sterilant / high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

Veltek Associates, Inc.
15 Lee Boulevard • Malvern, PA 19355-1234
Tel: (610)644-8335 Fax: (610)644-8336 E-mail: vai@sterile.com VL-501-98 12/18/03

DECON-CYCLE[®]

Sterile Pharmaceutical Clean Room Formula

Germicidal Detergent with Organic Soil Tolerance

MATERIAL SAFETY DATA SHEET

**COMPLIES WITH OSHA HAZARD COMMUNICATION STANDARD 29 CFR 1910.1200
(Complies with Commission Directive 91/155/EEC amended by 2001/58/EC)**



Veltek Associates, Inc.

15 Lee Boulevard • Malvern, PA 19355-1234

Tel: (610)644-8335 Fax: (610)644-8336 E-mail: vai@sterile.com VL-501-98 12/18/03

DECON-CYCLE[®]
Sterile Pharmaceutical Clean Room Formula
Germicidal Detergent with
Organic Soil Tolerance

MATERIAL SAFETY DATA SHEET
COMPLIES WITH OSHA HAZARD COMMUNICATION STANDARD 29 CFR
1910.1200
(Complies with Commission Directive 91/155/EEC amended by 2001/58/EC)



Veltek Associates, Inc.
15 Lee Boulevard • Malvern, PA 19355-1234
Tel: (610)644-8335 Fax: (610)644-8336 E-mail: vai@sterile.com VL-501-98 12/18/03

DECON-CYCLE®

MATERIAL SAFETY DATA SHEET

**COMPLIES WITH OSHA HAZARD COMMUNICATION STANDARD 29 CFR
1910.1200**

(Complies with Commission Directive 91/155/EEC amended by 2001/58/EC)

I. Identification of Substance:

Product Name: DECON-CYCLE®

MSDS #: DCY-98-01

Description: A concentrated low pH germicidal detergent effective in hard water and 5% organic soil load

MSDS Information	Manufacturer Information	Emergency Telephone Information
Date Prepared: November 6, 2002	Veltek Associates, Inc.	Notify your Supervisor
Prepared by:	15 Lee Boulevard	Emergencies Chemtrec (800) 424-9300
Art Vellutato, Jr.	Malvern, PA 19355	24-hour service
V.P. Technical Support Operations	TEL: (610)644-8335	
Reviewed by:		Effective Date: November 6, 2002
Art Vellutato, Sr.		Date Supersedes: 03/26/98
Technical Director		

II. Composition/Data on Components:

Hazardous Identity Components (Specific Chemical) Cas #'s OSHA PEL ACGIH TLV Other Limits % (Optional)

Hazardous Identity Components (Specific Chemical) Cas #'s OSHA PEL ACGIH TLV Other Limits % (Optional)

Isopropyl alcohol	67-63-0	400 ppm	400 ppm	N/E	
o-benzyl-p-chlorophenol	120-32-1	N/E	N/E	N/E	8.03%
2-Phenylphenol	90-43-7	N/E	N/E	N/E	6.28%
Dodecylbenzenesulfonic acid	21176-87-0	N/E	N/E	N/E	
Triethanolamine	102-21-6	N/E	N/E	N/E	
Hexylene Glycol	107-41-5	N/E	N/E	N/E	

Ingredients subject to reporting requirements of Section 313 of SARA Title III Ingredients listed as germicidal actives
N/A = Not Applicable N/E = Not Established N/D = Not Determined N/L = Not Listed

HIMIS RATING *Concentrate* HEALTH 2 REACTIVITY 0 FLAMMABILITY 2 PERSONAL B
HIMIS RATING *Diluted 1:256 (1/2 oz. per gal)* HEALTH 1 REACTIVITY 0 FLAMMABILITY 0 PERSONAL B

HAZARD RATINGS: 0 = Minimal 1 = Slight 2 = Moderate 3 = Serious 4 = Extreme A = Safety Glasses B = Safety Glasses & Gloves

Veltek Associates, Inc.

15 Lee Boulevard • Malvern, PA 19355-1234

Tel: (610)644-8335 Fax: (610)644-8336 E-mail: vai@sterile.com VL-501-98 12/18/03

O V E R V I E W
12

III. Hazards Identification:

HEALTH HAZARDS (Acute & Chronic)

ROUTES OF ENTRY:	Eyes?	X	Skin?	X	Ingestion?	X	Inhalation:	X
CARCINOGENICITY:	NTP?	NO	IARC Monographs?	NO	OSHA Regulated?			NO

HEALTH HAZARDS: Will cause eye irritation and skin irritation with excessive or prolonged exposure. Harmful or fatal if swallowed. Corrosive to lungs.

SIGNS AND SYMPTOMS OF EXPOSURE: If in eyes: Burning sensation, watering or redness. If on skin: Redness, irritation or burning sensation with prolonged exposure. If spray mist is inhaled: Lung irritation and/or burning sensation. If Swallowed: Drowsiness, irregular pulse, loss of consciousness.

MEDICAL CONDITIONS GENERALLY AGGRAVATED BY EXPOSURE: None known

IV. First Aid Measures:

EMERGENCY AND FIRST AID PROCEDURES:

If in eyes: Flush with water for 15 minutes.

If on skin: Wash liberally with soap and water.

If on clothes: Wash before reuse.

If swallowed: Dilute by drinking milk, egg whites, or large quantities of water. Do not induce vomiting without medical advice.

If spray mist is inhaled: Remove to fresh air.

Always get medical attention when product is swallowed or when symptoms are significant or persist

V. Firefighting Measures:

FLASH POINT 107^o F (Combustible Liquid) (METHOD USED): SETA Flammable Limits: N/E LEL: N/E UEL: N/E

EXTINGUISHING MEDIA: In case of fire use water spray, CO₂, Dry Chemical or Foam

SPECIAL FIRE FIGHTING PROCEDURES: Use water spray to cool fire-exposed containers and to protect personnel

UNUSUAL FIRE AND EXPLOSION HAZARDS: Burning may release combustion products such as Phenolics, carbon monoxide, carbon dioxide, and chlorine

VI. Accidental Release Measures:

STEPS TO BE TAKEN IN CASE MATERIAL IS RELEASED OR SPILLED: Small spills (less than 1 gallon) may be washed down a drain or cleaned up and disposed of into a sanitary sewer system. Large spills (more than 1 gallon) should be contained and collected (by absorption or vacuuming) then disposed of properly.

VII. Handling and Storage:

STORAGE & HANDLING PRECAUTIONS: KEEP OUT OF REACH OF CHILDREN. Do not contaminate water, food, or feed. Do not reuse empty containers. Keep containers clean and closed.

Veltek Associates, Inc.

15 Lee Boulevard • Malvern, PA 19355-1234

Tel: (610)644-8335 Fax: (610)644-8336 E-mail: vai@sterile.com VL-501-98 12/18/03

O V E R V I E W
13

VIII. Exposure Control and Personal Protection:

RESPIRATORY PROTECTION: General ventilation is normally adequate. Protect breathing system with an appropriate respirator if spraying causes exposure to mists.

GLOVES: Required for exposure to the concentrate when diluting or for long exposures to end-use dilutions. Persons sensitive to cleaning chemicals should always wear gloves.

EYE PROTECTION: Use Safety Glasses when diluting concentrates or when splashing or spraying of diluted product into the eyes is likely.

OTHER PROTECTIVE CLOTHING OR EQUIPMENT: Not normally necessary for routine handling.

WORK/HYGIENE PRACTICES: Wash hands after using. Do not get into eyes, on skin, or clothing. May be harmful if swallowed. Protect food and drink from contamination by product.

VENTILATION: Local Exhaust: Normally Sufficient Mechanical (General): Not Normally Needed

NOTE: *The information herein is presented in good faith and believed to be accurate as of the effective date shown above. However, no warranty, expressed or implied is given. It is the responsibility of the user to comply with all Federal, State and Local laws and regulations.*

IX. Physical and Chemical Properties:

Appearance and Odor: Clear Amber Liquid (that may darken with age) with a slight alcohol odor

Boiling Point:	212 ^o F	Specific Gravity: @25 ^o C	1.092 +/- .01	pH (Concentrate):	2.00 +/- 0.50
Vapor Pressure	N.E.	Solubility in Water:	Complete	pH (Diluted): 1:256:	3.75 +/- 0.75
Vapor Density:	N.E.	Melting Point:	N/A	Evaporation Rate:	N.E.

X. Stability and Reactivity:

STABILITY: Stable

HAZARDOUS POLYMERIZATION: Will Not Occur

INCOMPATIBILITY (Materials or Conditions to Avoid): Extreme high temperatures or open flame

HAZARDOUS DECOMPOSITION OR BYPRODUCTS: None known

XI. Toxicological Information:

LD ₅₀ Oral Toxicity in Rats:	LD ₅₀ >500 mg/kg and <5000 mg/kg
LD ₅₀ Dermal Toxicity in Rabbits:	LD ₅₀ > 2.0 g/kg of body weight
Primary Dermal Irritation in Rabbits:	Corrosive (Primary Irritation Index: 4.33)
Primary Eye Irritation/corrosion in rabbits:	Irritant
Inhalation Toxicity in Rats:	LC ₅₀ > 2.3 mg/L
Delayed Contact Dermal Sensitization:	Not a Sensitizer

XII. Ecological Information:

EPA Reg. 211-62 is Readily Biodegradable: 74% within 28 days.

Consult local authorities pertaining to the release of material to the environment.

XIII. Disposal Considerations:

WASTE DISPOSAL METHOD: Dispose of product by using in accordance with label directions or in accordance with local, state, and federal regulations. This product is NOT considered a RCRA Hazardous Waste as identified in Federal Regulations 40 CFR 261.

XIV. Transport Information:

Flammable Liquid, Corrosive N.O.S. (Isopropanol, Citric Acid)

Class 3 (Flammable), Packaging Group III, UN 2924

Veltek Associates, Inc.

15 Lee Boulevard • Malvern, PA 19355-1234

Tel: (610)644-8335 Fax: (610)644-8336 E-mail: vai@sterile.com VL-501-98 12/18/03

O V E R V I E W

14

XV. Regulations:

PLEASE CONSULT YOUR STATE AND LOCAL AUTHORITY FOR INGREDIENTS LISTED AS PERTAINING TO REPORTING REQUIREMENTS FOR HAZARDOUS SUBSTANCE LIST.

No special regulation other than transportation and disposal (see section VII and section IX). Conform to Federal, State and Local Regulations for disposal. Follow usage instructions on container.

XVI. Other Information:

NOTE: The information herein is presented in good faith and believed to be accurate as of the effective date shown above. However, no warranty, expressed or implied is given. It is the responsibility of the user to comply with all Federal, State and Local laws and regulations.

Veltek Associates, Inc.

15 Lee Boulevard • Malvern, PA 19355-1234

Tel: (610)644-8335 Fax: (610)644-8336 E-mail: vai@sterile.com VL-501-98 12/18/03

DECON-CYCLE[®]

**Sterile Pharmaceutical Clean Room Formula
Low pH Phenolic Germicidal Detergent**

Lot Specific Sterile Documentation

(received with each shipment)

**Certificate of Analysis
Certificate of Sterility
Certificate of Irradiation**

(Please contact VAI for a sample of this documentation)