



BioClean Extra™

EXTRA LENGTH STERILE LATEX CLEANROOM GLOVES

BLAS



Bioclean Extra™ sterile 400mm (16") elbow length latex cleanroom gloves have been developed for users who require increased protection and control against product contamination.

They are textured for excellent grip, particularly when wet, and have extremely low levels of particulates and extractables.

Sterile cleanroom gloves are packed in 'Easy-tear'™ plastic packaging that is silicone, amine and amide-free to ensure compatibility with ISO Class 4 & EU GMP Grade A environments.



Low levels of latex proteins minimise risk of potential allergies and the extra length provides improved user protection. These latex cleanroom gloves are also available in a non-sterile version and in a range of sizes to suit all users. They also comply to the European PPE Category 3 Complex Design.

"Extra Length Protection!"

NITRITEX LTD.	T: +44 (0) 1638 663338	E: info@nitritex.com
NITRITEX ASIA PTE LTD.	T: +65 6536 9123	E: nitritex@pacific.net.sg
NITRITEX PUERTO RICO Inc.	T: +1 787 232 5084	E: info@nitritex.com.pr
NITRITEX CANADA LTD.	T: +1 (905) 946 9539	E: info@nitritexcanada.com
NITRITEX MALAYSIA Sdn. Bhd.	T: +(603) 3122 2614	E: operations@nitritex.com.my





BioClean Extra™

EXTRA LENGTH STERILE LATEX CLEANROOM GLOVES

BLAS

Colour:	Natural.
Length:	400mm (16").
Material:	Natural rubber latex.
Type:	Ambidextrous.
Surface:	Textured.
Other Features:	Powderfree, beaded cuff.
Protein Level:	Less than 50µg/gram.
Sterilisation Method:	Gamma Irradiation.
Packing:	Bagged.
AQL:	1.5 for freedom from pinholes (1000ml) water test.
Particulates:	Processed in a cleanroom environment to ensure compatibility with ISO Class 4 operating conditions.
Extractables:	Extended leaching in pure de-ionised water ensures low levels of chlorides and extractable metals.
Physical Properties:	Manufactured to exceed the requirements of the PPE Gloves Standard EN420:2003 and the European Medical Gloves Standard EN455 (Parts 1, 2 & 3).
Quality Systems:	Manufactured in a facility holding ISO 9001:2000 and 89/686/EEC Article 11b certifications.
Sampling & Inspection:	In accordance with ISO 2859.
Storage Instructions:	Store in a cool, dry place (5°C - 35°C) away from direct sunlight and fluorescent light.
Shelf Life:	Five (5) years from date of manufacture.
Bar Codes:	Each size carries a unique bar code on the carton (ITF-14). This allows automatic data entry into Stock Control Systems.
Lot Number:	This is shown on the sealed pouch and outer carton to provide product traceability.
Packaging:	One cuffed pair per inner plastic wallet; One inner wallet per 'Easy Tear' PE bag; 10 pairs per sealed outer PE bag; 20 outer bags per carton liner; One carton liner per carton (200 pairs).
Compatibility:	ISO Class 4 & EU GMP Grade A.
Regulatory Classification:	Category 3 - Complex Design Personal Protective Equipment conforming to EN420:2003 & EN374-1 & 2:2003.
Product Codes and Sizes:	BLAS-S: Small BLAS-M: Medium BLAS-L: Large BLAS-XL: Extra Large

"Extra Length Protection!"

BLAS/p/DS7



GB04/61092

