



ΩMEGA Pro-tech XP™

EXTRA PROTECTION NON-STERILE NITRILE GLOVES

OPT XP



Omega pro-tech XP™ Non-sterile powderfree ambidextrous nitrile gloves deliver excellent protection against bloodborne pathogens and are 100% Latex free, thereby eliminating all risk of Type 1 latex allergic reactions.

Non-sterile powderfree ambidextrous nitrile gloves are 300mm (12") in length to provide extra protection in the cuff area. These violet coloured nitrile gloves are resistant to a wide range of chemicals, fatty acids and bodily fluids. They have strong puncture resistance with no compromise on tactile performance and dexterity.

Our ambidextrous nitrile gloves have a beaded cuff which ensures easy donning and the textured surface provides an excellent grip. The low modulus nitrile compound quickly and comfortably conforms to the hand giving excellent tactility.

These gloves are available in a wide range of sizes to suit all users and the perforated dispenser box provides easy dispensing. These gloves comply to the European PPE Category 3 Complex Design.



"Extra Protection & Chemical Resistance!"

NITRITEX LTD.	T: +44 (0) 1638 663338	E: info@nitritex.com
NITRITEX ASIA PTE LTD.	T: +65 6536 9123	E: nitritex@pacific.net.sg
NITRITEX PUERTO RICO Inc.	T: +1 787 232 5084	E: info@nitritex.com.pr
NITRITEX CANADA LTD.	T: +1 (905) 946 9539	E: info@nitritexcanada.com
NITRITEX MALAYSIA Sdn. Bhd.	T: +(603) 3122 2614	E: operations@nitritex.com.my





ΩMEGA Pro-tech XP™

EXTRA PROTECTION NON-STERILE NITRILE GLOVES

OPT XP

Colour:	Violet.
Length:	300mm (12").
Material:	Nitrile.
Type:	Ambidextrous.
Surface:	Finger textured.
Other Features:	Powderfree, latex free, beaded cuff.
Protein Level:	Zero. Nitrile contains no natural latex proteins.
Packing:	Boxed.
AQL:	1.5 for holes (1000ml water test).
Physical Properties:	Manufactured to exceed the requirements of the PPE Gloves Standard EN420:2003 and the European Medical Gloves Standard EN455 (Parts 1, 2 & 3).
Quality Systems:	Manufactured in a facility holding ISO 9001:2000 and 89/686/EEC Article 11b certifications.
Sampling & Inspection:	In accordance with ISO 2859.
Storage Instructions:	Store in a cool, dry place (5°C - 35°C) away from direct sunlight and fluorescent light.
Shelf Life:	Five (5) years from date of manufacture.
Bar Codes:	Each size carries a unique bar code on the inner dispenser box (EAN-13) and the outer carton (ITF-14). This allows automatic data entry into Stock Control Systems.
Lot Number:	This appears on the inner dispenser box and the outer carton to provide traceability.
Packaging:	100 pieces per inner dispenser box; 10 dispenser boxes per carton (1000 pieces).
Regulatory Classification:	Category 3 - Complex Design Personal Protective Equipment conforming to EN420:2003 & EN374-1 & 2:2003.
Product Codes and Sizes:	OPT XP-XS: X-Small OPT XP-S: Small OPT XP-M: Medium OPT XP-L: Large OPT XP-XL: X-Large OPT XP-XXL: XX-Large

"Extra Protection & Chemical Resistance!"

OPTXP/PS IN



GB04/61092

