



BNPS

BioClean N-Plus™

STERILE 400mm (16") NITRILE CLEANROOM GLOVES



BioClean N-Plus™ Sterile 400mm (16") Length Nitrile Cleanroom Gloves have been specifically developed to provide extra protection above the wrist.

Formulated for superior chemical resistance **BioClean N-Plus™** Sterile 400mm (16") Length Powder-free Nitrile Cleanroom Gloves are compatible for use in ISO Class 5 cleanrooms. They are hand-specific and anatomically shaped to give excellent user comfort.

BioClean N-Plus™ Sterile 400mm (16") Length Nitrile Cleanroom Gloves have a textured surface to ensure good grip, even in wet conditions. These gloves are available in a wide range of sizes to suit all users.

"Full Nitrile Protection!"

NITRITEX ASIA PTE LTD
11 Keng Cheow Street
#03-09 Riverside Piazza
Singapore 059608
Republic of Singapore
Tel: +65 6536 9123
Fax: +65 6438 5159
Email: nitritex@pacific.net.sg

NITRITEX CANADA LTD
7030 Woodbine Avenue
Suite 500, Markham
Ontario L3R 6G2
Canada
Tel: +1 (905) 946 9539
Fax: +1 (905) 946 8584
Email: info@nitritexcanada.com

NITRITEX LTD
Minton Enterprise Park
Oaks Drive, Newmarket
Suffolk, CB8 7YY
United Kingdom
Tel: +44 (0) 1638 663338
Fax: +44 (0) 1638 668890
Email: sales@nitritex.com





BioClean N-PlusTM

STERILE 400mm (16") NITRILE CLEANROOM GLOVES

- Manufacture:** Nitritex Ltd.
- Features:** Hand-specific, beaded cuff, powder-free.
- Surface:** Textured.
- Colour:** White.
- Length:** 400mm.
- Material:** Nitrile.
- Protein Level:** Zero. Nitrile contains no natural latex proteins.
- Physical Properties:** **BioClean N-PlusTM** Sterile Cleanroom Gloves are manufactured to exceed the requirements of the European Medical Gloves Standard EN455 (Parts 1 & 2).
- Quality Systems:** Manufactured in a facility holding the following registration: ISO 9001:2000
- Sampling & Inspection:** In accordance with ISO 2859.
- AQL:** 1.5 for holes (1000ml water test).
- Storage Instructions:** Store in a cool dry place away from direct sunlight and heat (5°C-35°C).
- Lot Number:** This appears on both the sealed pair pouch and the outer carton to provide traceability.
- Bar Codes:** Each size carries a unique bar code on the outer carton (ITF 14) to allow automatic entry into Stock Control Systems.
- Packaging:** 1 pair per inner wallet/Easy-tear bag;
10 pairs per sealed inner bag;
20 inner bags per carton liner;
One carton liner per carton (200 pairs).
- Product Codes:** Each size is given an individual code for re-order purposes as follows:

SIZE	Re-order Code
6.0	BNPS 60
6.5	BNPS 65
7.0	BNPS 70
7.5	BNPS 75
8.0	BNPS 80
8.5	BNPS 85
9.0	BNPS 90
10.0	BNPS 100

"Full Nitrile Protection!"



GB04/61092