



BioClean Alpha™

STERILE POWDERFREE LATEX CLEANROOM GLOVES

AL300T



BioClean Alpha™ (AL300T) sterile 300mm (12") length latex powderfree cleanroom gloves have been developed for use as the '*outer pair*' of gloves in combination with **BioClean Emerald™** or **BioClean Jade™** gloves or as a standard pair of gloves when '*double gloving*' is not required. The extra strength cuff prevents tearing whilst donning.

These sterile 300mm (12") length latex powderfree cleanroom gloves have a textured surface which provides an excellent grip.

Our cleanroom gloves are hand specific, anatomically shaped for comfort and are available in a wide range of sizes to suit all users.

This is a specially developed version for those users requiring greater tack.

"Extra Strength, Tack & Quality!"

NITRITEX LTD.	T: +44 (0) 1638 663338	E: info@nitritex.com
NITRITEX ASIA PTE LTD.	T: +65 6536 9123	E: nitritex@pacific.net.sg
NITRITEX PUERTO RICO Inc.	T: +1 787 232 5084	E: info@nitritex.com.pr
NITRITEX CANADA LTD.	T: +1 (905) 946 9539	E: info@nitritexcanada.com
NITRITEX MALAYSIA Sdn. Bhd.	T: +(603) 3122 2614	E: operations@nitritex.com.my


www.nitritex.com



BioClean Alpha™

STERILE POWDERFREE LATEX CLEANROOM GLOVES

AL300T

Colour:	Natural.
Length:	300mm (12").
Material:	Natural rubber latex.
Type:	Hand specific.
Surface:	Textured.
Other Features:	Powderfree, beaded cuff, anatomical shape.
Protein Level:	Less than 50µg/gram.
Sterilisation Method:	Gamma Irradiation.
Packing:	Bagged.
AQL:	1.5 for freedom from pinholes (1000ml) water test.
Particulates:	Processed in a cleanroom environment to ensure compatibility with ISO Class 4 operating conditions.
Extractables:	Extended leaching in pure de-ionised water ensures low levels of chlorides and extractable metals.
Physical Properties:	Manufactured to exceed the requirements of the PPE Gloves Standard EN420:2003 and the European Medical Gloves Standard EN455 (Parts 1, 2 & 3).
Quality Systems:	Manufactured in a facility holding ISO 9001:2000 and 89/686/EEC Article 11b certifications.
Sampling & Inspection:	In accordance with ISO 2859.
Storage Instructions:	Store in a cool, dry place (5°C - 35°C) away from direct sunlight and fluorescent light.
Shelf Life:	Five (5) years from date of manufacture.
Bar Codes:	Each size carries a unique bar code on the carton (ITF-14). This allows automatic data entry into Stock Control Systems.
Lot Number:	This is shown on the sealed pouch and outer carton to provide product traceability.
Packaging:	One cuffed pair per inner plastic wallet; One inner wallet per 'Easy Tear' PE bag; 10 pairs per sealed outer PE bag; 20 outer bags per carton liner; One carton liner per carton (200 pairs).
Compatibility:	ISO Class 4 & EU GMP Grade A.
Regulatory Classification:	Category 3 - Complex Design Personal Protective Equipment conforming to EN420:2003 & EN374-1 & 2:2003.
Product Codes and Sizes:	AL300T-60: 6.0 AL300T-65: 6.5 AL300T-70: 7.0 AL300T-75: 7.5 AL300T-80: 8.0 AL300T-85: 8.5 AL300T-90: 9.0 AL300T-100: 10.0

"Extra Strength, Tack & Quality!"

AL300/PDS6



GB04/61092

