



BioClean Maxima™

BLLS

EXTRA LENGTH STERILE LATEX CLEANROOM GLOVES



BioClean Maxima™

Extra Length Sterile 600mm

(24") Length Latex Powder-free Cleanroom Gloves have been developed to provide the ultimate in user protection.

BioClean Maxima™ Extra Length Sterile 600mm (24") Length Latex Powder-free Cleanroom Gloves have excellent grip and are particularly suitable for use in environments where fluid contamination is a risk. They are compatible for use in ISO Class 5 cleanrooms.

BioClean Maxima™ Extra Length Sterile 600mm (24") Length Latex Powder-free Cleanroom Gloves are hand specific and available in a wide range of sizes to suit all users.

"The Ultimate in User Protection"

NITRITEX ASIA PTE LTD
11 Keng Cheow Street
#03-09 Riverside Piazza
Singapore 059608
Republic of Singapore
Tel: +65 6536 9123
Fax: +65 6438 5159
Email: nitritex@pacific.net.sg

NITRITEX CANADA LTD
7030 Woodbine Avenue
Suite 500, Markham
Ontario L3R 6G2
Canada
Tel: +1 (905) 946 9539
Fax: +1 (905) 946 8584
Email: info@nitritexcanada.com

NITRITEX LTD
Minton Enterprise Park
Oaks Drive, Newmarket
Suffolk, CB8 7YY
United Kingdom
Tel: +44 (0) 1638 663338
Fax: +44 (0) 1638 668890
Email: sales@nitritex.com





BioClean Maxima™

EXTRA LENGTH STERILE LATEX CLEANROOM GLOVES

- Manufacturer:** Nitritex Ltd
- Features:** Powderfree, beaded cuff, hand specific, anatomical shape
- Surface:** Textured
- Colour:** White
- Length:** 600mm (24")
- Material:** Natural rubber latex
- Protein Level:** Less than 50 µg/gram
- Sterilisation Method:** Gamma Irradiation
- Physical Properties:** **BioClean Maxima™** sterile latex cleanroom gloves are manufactured to exceed the requirements of the European Medical Gloves Standard (EN455 Parts 1,2 & 3)
- Particulates:** **BioClean Maxima™** sterile latex cleanroom gloves are processed in a cleanroom environment to ensure compatibility with ISO Class 5 operating conditions.
- Extractables:** Extended leaching in pure de-ionised water ensures extremely low levels of chlorides and extractable metals.
- Packaging:** One cuffed pair per inner plastic wallet;
One inner wallet per 'Easy Tear' PE bag;
25 pairs per inner box;
4 inner boxes per carton (100 pairs)
- Quality Systems:** Manufactured in a facility holding the following registration:
ISO 9001:2000
- Sampling & Inspection:** In accordance with ISO 2859
- AQL:** 1.5 for freedom from pinholes (1000ml water test)
- Storage Instructions:** Store in a cool (5°C - 35°C), dry place away from direct sunlight and fluorescent light.
- Shelf Life:** Five (5) years from date of manufacture
- Lot Number:** This is shown on the sealed 'Easy Tear' bag, the inner box and the carton to provide product traceability
- Bar Codes:** Each size carries a unique bar code on the carton (ITF-14).
This allows automatic data entry into stock control systems

Product Codes:	Size	Code
	6.0	BLLS-60
	6.5	BLLS-65
	7.0	BLLS-70
	7.5	BLLS-75
	8.0	BLLS-80
	8.5	BLLS-85
	9.0	BLLS-90

"The Ultimate in User Protection"

