




BUPS

# BioClean Ultimate™

## STERILE 300mm (12") LENGTH POWDERFREE SYNTHETIC CLEANROOM GLOVES



**BioClean Ultimate™** Sterile Synthetic Cleanroom Gloves are manufactured from a proprietary co-polymer blend. With no latex proteins Type 1 latex allergic reactions are eliminated.

With their anatomical hand-specific shape, **BioClean Ultimate™** Sterile Synthetic Cleanroom Gloves have been designed for those looking for the ultimate in tactility. They have the fit and feel of natural latex gloves and the beaded cuff ensures a secure fit at the wrist. The textured surface provides a good grip even in wet environments.

Available in half sizes from 6.0 to 9.0 to suit all users, **BioClean Ultimate™** Sterile Synthetic Cleanroom Gloves are 300mm (12") in length for maximum protection in the cuff area. They are gamma sterilised and are compatible for use in sterile ISO Class 5 cleanroom environments.

**"The Ultimate Glove!"**

**NITRITEX ASIA PTE LTD**  
11 Keng Cheow Street  
#03-09 Riverside Piazza  
Singapore 059608  
Republic of Singapore  
Tel: +65 6536 9123  
Fax: +65 6438 5159  
Email: nitritex@pacific.net.sg

**NITRITEX CANADA LTD**  
7030 Woodbine Avenue  
Suite 500, Markham  
Ontario L3R 6G2  
Canada  
Tel: +1 (905) 946 9539  
Fax: +1 (905) 946 8584  
Email: info@nitritexcanada.com

**NITRITEX LTD**  
Minton Enterprise Park  
Oaks Drive, Newmarket  
Suffolk, CB8 7YY  
United Kingdom  
Tel: +44 (0) 1638 663338  
Fax: +44 (0) 1638 668890  
Email: sales@nitritex.com





# BioClean Ultimate™

## STERILE 300mm (12") LENGTH POWDERFREE SYNTHETIC CLEANROOM GLOVES

- Manufacturer:** Nitritex Ltd.
- Features:** Hand-specific, beaded cuff, powderfree.
- Surface:** Textured.
- Colour:** Natural.
- Length:** 300mm (12").
- Material:** Proprietary Co-polymer blend
- Protein Level:** Contains no natural latex proteins.
- Physical Properties:** **BioClean Ultimate™** cleanroom gloves are manufactured to exceed the requirements of the European Medical Gloves Standard EN455 (Parts 1 & 2).
- Quality Systems:** Manufactured in a facility holding the following registration: ISO 9001:2000.
- Sampling & Inspection:** In accordance with ISO 2859.
- AQL:** 1.5 for holes (1000ml water test).
- Storage Instructions:** Store in a cool dry place away from direct sunlight and heat (5°C-35°C).
- Lot Number:** This appears on both the sealed pair pouch and the outer carton to provide traceability.
- Bar Codes:** Each size carries a unique bar code on the outer carton (ITF 14) to allow automatic entry into Stock Control Systems.
- Packaging:** One cuffed pair per inner PE wallet;  
One inner wallet per 'Easy-Tear' PE pouch;  
10 pairs per sealed inner bag;  
20 inner bags per carton liner;  
One carton liner per carton (200 pairs).
- Product Codes:** Each size is given an individual code for re-order purposes as follows:

SIZE	RE-ORDER CODE	SIZE	RE-ORDER CODE
6.0	BUPS 60	8.0	BUPS 80
6.5	BUPS 65	8.5	BUPS 85
7.0	BUPS 70	9.0	BUPS 90
7.5	BUPS 75		



GB04/61092



### "The Ultimate Glove!"