



BNAL

BioClean 100™

AMBIDEXTROUS EXTRA PROTECTION 400mm (16") LENGTH NITRILE CLEANROOM GLOVES



BioClean 100™ Non-sterile, Powder-free, Ambidextrous, Extra protection 400mm (16") Length Nitrile Cleanroom Gloves set new standards of excellence in cleanroom gloves. Double bagged for critical environment use, they are thoroughly washed in de-ionised water for ISO Class 5 compatibility.

BioClean 100™ Extra Protection Nitrile Cleanroom Gloves are 100% latex free, eliminating all latex allergy problems. They are fingertip textured for improved grip and for the safe handling of products, even when wet.

They are resistant to a wide range of chemicals, stronger than natural latex gloves and have excellent ESD properties. The 400mm (16") length provides added user protection above the cuff area.

BioClean 100™ Non-sterile, Powder-free, Ambidextrous, Extra protection Nitrile Cleanroom Gloves are available in a wide range of sizes to suit all users.

"Latex Free Extra Protection!"

NITRITEX ASIA PTE LTD
11 Keng Cheow Street
#03-09 Riverside Piazza
Singapore 059608
Republic of Singapore
Tel: +65 6536 9123
Fax: +65 6438 5159
Email: nitritex@pacific.net.sg

NITRITEX CANADA LTD
7030 Woodbine Avenue
Suite 500, Markham
Ontario L3R 6G2
Canada
Tel: +1 (905) 946 9539
Fax: +1 (905) 946 8584
Email: info@nitritexcanada.com

NITRITEX LTD
Minton Enterprise Park
Oaks Drive, Newmarket,
Suffolk, CB8 7YY
United Kingdom
Tel: +44 (0) 1638 663338
Fax: +44 (0) 1638 668890
Email: sales@nitritex.com





BioClean 100™

AMBIDEXTROUS EXTRA PROTECTION 400mm (16") LENGTH NITRILE CLEANROOM GLOVES

Manufacture:	Nitritex Ltd
Features:	Powder-free, ambidextrous, beaded cuff
Surface:	Fingertip textured
Material:	Nitrile
Physical Properties:	Manufactured to exceed the requirements of European Medical Gloves Standard EN 455 (Parts 1 & 2)
Colour:	White
Length:	400mm (16")
Particulates:	Processed in a cleanroom environment to ensure compatibility with ISO Class 5 operating conditions
Protein Level:	Zero. Nitrile contains no natural latex proteins
AQL:	1.5 for holes (1000ml water test)
Extractables:	Extended leaching in de-ionised water ensures extremely low levels of chlorides and extractable metals
Registering Authority:	SGS UK Limited
Quality Systems:	Manufactured in a facility holding the ISO 9001:2000 registration
Sampling & Inspection:	In accordance with ISO 2859
Storage Instructions:	Store in a cool dry place (5°C-35°C) away from direct sunlight and heat
Lot Number:	This appears on the sealed polybag, and the outer carton to provide traceability
Packing:	100 pieces per sealed inner polybag; one inner polybag per sealed outer polybag; 10 outer polybags per carton (1000 pieces)

Product Codes:	SIZE	RE-ORDER CODE
	Small	BNAL-S
	Medium	BNAL-M
	Large	BNAL-L
	X-Large	BNAL-XL

"Latex Free Extra Protection!"

