



BioClean Jade™

STERILE NITRILE GLOVES

BJAS



BioClean Jade™ Sterile Powderfree Ambidextrous Nitrile Gloves are one of the most technically advanced gloves available today.

They are 100% sulphur and accelerator-free, thereby virtually eliminating Type 4 allergies. Whilst their smooth surface makes them ideal as an underglove where double-donning is required, they can be used by themselves when extra sensitivity is required. Compatible for use in sterile ISO Class 4 and EU GMP Grade A environments, they are 100% latex free, eliminating Type 1 latex allergies.

With a super-soft formulation, **BioClean Jade™** Sterile Nitrile Gloves are 300mm (12") in length and are available in a wide range of sizes to fit all users.

"The Technically Advanced Underglove!"

NITRITEX LTD. Tel: +44 (0) 1638 663338 E-mail: info@nitritex.com

NITRITEX ASIA PTE LTD. Tel: +65 6536 9123 E-mail: nitritex@pacific.net.sg

NITRITEX CANADA LTD. Tel: +1 (905) 946 9539 E-mail: info@nitritexcanada.com

NITRITEX MALAYSIA Sdn. Bhd. Tel: (603) 3122 2614 E-mail: operations@nitritex.com.my

 **nitritex**
www.nitritex.com



BioClean Jade™

STERILE NITRILE GLOVES

BJAS

Manufacture:	Nitritex Limited
Features:	Beaded cuff, powderfree, ambidextrous, smooth surface
Colour:	Green
Length:	300mm (12")
Material:	Nitrile
Protein Level:	Zero. Nitrile contains no natural latex proteins
Accelerators:	Zero. 100% accelerator-free
Physical Properties:	BioClean Jade™ Nitrile cleanroom gloves are manufactured to exceed the requirements of the European Medical Gloves Standard EN455 (Parts 1 & 2)
Quality Systems:	Manufactured in a facility holding the ISO 9001:2000 registration
Sampling & Inspection:	In accordance with ISO 2859
AQL:	1.5 for holes (1000ml water test)
Sterilization Method:	Gamma Irradiation
Storage Instructions:	Store in a cool dry place away from direct sunlight and heat (5°C-35°C)
Shelf Life:	Five years from date of manufacture
Lot Number:	This appears on the pair pouch and carton to provide traceability
Bar Codes:	Each size carries a unique bar code on the carton (ITF 14) to allow automatic entry into Stock Control Systems
Packaging:	One cuffed pair per inner plastic wallet; One inner wallet per sealed Easy-tear bag; 10 Easy-tear bags per outer PE bag; 20 outer PE bags per carton (200 pairs)
Product Codes:	Each size is given an individual code for re-order purposes as follows:

Size	Re-order Code
5.0-5.5	BJAS 5055
6.0-6.5	BJAS 6065
7.0-7.5	BJAS 7075
8.0-8.5	BJAS 8085
9.0	BJAS 9090
10.0	BJAS 1010

"The Technically Advanced Underglove!"

BJAS/PDS 1



GB04/61092

