




NAN4

BioClean Nano 4™

NON-STERILE NITRILE AMBIDEXTROUS CLEANROOM GLOVES



BioClean Nano 4™

The BioClean Nano 4™ Non-sterile, powderfree, ambidextrous glove from Nitritex, sets a new standard of excellence in Nitrile cleanroom gloves.

Design

The BioClean Nano 4™ with lightly textured fingers for improved grip, have an excellent tack for the safe handling of products, even when wet. They are resistant to a wide range of chemicals, have a far greater puncture resistance than natural latex gloves and have excellent ESD properties. The 300mm (12") length provides added user protection in the cuff area.

Allergens

BioClean Nano 4™ gloves are 100% latex free, eliminating latex allergy problems.

Packaging

The gloves are double bagged for critical environment use and are thoroughly washed in pure de-ionised water for ISO Class 4 compatibility.

This product is available in a wide range of sizes to suit all users.

"A New Standard of Excellence!"

NITRITEX LTD. Minton Enterprise Park, Oaks Drive, Newmarket, Suffolk, CB8 7YY, United Kingdom
Tel: +44 (0) 1638 663338 Fax: +44 (0) 1638 668890 E-mail: sales@nitritex.com

NITRITEX CANADA LTD. 7030 Woodbine Avenue, Suite 500, Markham, Ontario, L3R 6G2, Canada
Tel: +1 (905) 946 9539 Fax: +1 (905) 946 8584 E-mail: info@nitritexcanada.com

NITRITEX ASIA PTE LTD. 11 Keng Cheow Street, #03-09 Riverside Piazza, Singapore 059608, Republic of Singapore
Tel: +65 6536 9123 Fax: +65 6438 5159 E-mail: nitritex@pacific.net.sg

 **nitritex**
www.nitritex.com



BioClean Nano 4™

NAN4

NON-STERILE NITRILE AMBIDEXTROUS CLEANROOM GLOVES

Manufacture:	Nitritex Ltd														
Features:	Powderfree, ambidextrous shape, beaded cuff														
Surface:	Textured fingers														
Colour:	White														
Length:	300mm (12")														
Material:	Nitrile (synthetic latex)														
Protein Level:	Zero. Nitrile contains no natural latex proteins														
Physical Properties:	BioClean Nano 4™ non-sterile powderfree ambidextrous nitrile cleanroom gloves are manufactured to exceed the requirements of the European Medical Glove Standard (EN455 Parts 1 & 2)														
Particulates:	BioClean Nano 4™ gloves are processed in a cleanroom environment to ensure compatibility with ISO Class 4 operating conditions														
Extractables:	Extended leaching in pure de-ionised water ensures extremely low levels of chlorides and extractable metals														
Packaging:	100 pieces per sealed inner polybag; one inner polybag per sealed outer polybag; 10 outer polybags per carton (1000 pieces)														
Quality Systems:	Manufactured in a facility holding the following registration: ISO 9001:2000														
Sampling & Inspection:	In accordance with ISO 2859														
AQL:	1.5 for holes (1000ml water test)														
Storage Instructions:	Store in a cool dry place (5°C - 35°C) away from heat and direct sunlight														
Lot Number:	This appears on the sealed outer bag and the carton to provide traceability														
Bar Codes:	Each size carries a unique bar code on the outer bag of 100 pieces (EAN-13) and the outer carton (ITF-14). This allows automatic data entry into Stock Control Systems														
Product Codes:	<table border="1"> <thead> <tr> <th>Size</th> <th>Code</th> </tr> </thead> <tbody> <tr> <td>X-Small</td> <td>NAN4-XS</td> </tr> <tr> <td>Small</td> <td>NAN4-S</td> </tr> <tr> <td>Medium</td> <td>NAN4-M</td> </tr> <tr> <td>Large</td> <td>NAN4-L</td> </tr> <tr> <td>X-Large</td> <td>NAN4-XL</td> </tr> <tr> <td>XX-Large</td> <td>NAN4-XXL</td> </tr> </tbody> </table>	Size	Code	X-Small	NAN4-XS	Small	NAN4-S	Medium	NAN4-M	Large	NAN4-L	X-Large	NAN4-XL	XX-Large	NAN4-XXL
Size	Code														
X-Small	NAN4-XS														
Small	NAN4-S														
Medium	NAN4-M														
Large	NAN4-L														
X-Large	NAN4-XL														
XX-Large	NAN4-XXL														

"A New Standard of Excellence!"

Design & artwork by John Tunon/KCSD
Graphics & Media Manager
- Nitritex Ltd -
NAN4/PD51



GB04/61092

