



BDBS-G

BioClean DB™

STERILE POUCH STYLE MASK WITH NECK GUARD



BioClean DB™ Sterile Pouch Style Mask with neck guard comes from the renowned **BioClean** range of products.

The **BioClean DB™** Sterile Pouch Style Mask with neck guard has an innovative, patented spacer which maintains a clear breathing space for the wearer, improving the overall comfort and ease of breathing. It also has an integral neck guard to close the breach in protection between the coverall and throat of the wearer.

The **BioClean DB™** Sterile Pouch Style Mask with neck guard has knitted headbands and all seams are ultrasonically bonded. The plastic coated, soft malleable steel noseband is fully enclosed.

The **BioClean DB™** Sterile Pouch Style Mask with neck guard inner and outer facing layers are constructed of spunbond non-woven polyester with the middle filter layer of meltblown polyester. The spacer is a non-latex band and the neck guard is spunbond non-woven polyester.

The **BioClean DB™** Sterile Pouch Style Mask with neck guard is supplied Gamma Sterilised and is suitable for use in ISO Class 4 and higher environments.

"Comfort and Protection!"

NITRITEX ASIA PTE LTD
11 Keng Cheow Street
#03-09 Riverside Piazza
Singapore 059608
Republic of Singapore
Tel: +65 6536 9123
Fax: +65 6438 5159
Email: nitritex@pacific.net.sg

NITRITEX CANADA LTD
7030 Woodbine Avenue
Suite 500, Markham
Ontario L3R 6G2
Canada
Tel: +1 (905) 946 9539
Fax: +1 (905) 946 8584
Email: info@nitritexcanada.com

NITRITEX LTD
Minton Enterprise Park
Oaks Drive, Newmarket
Suffolk, CB8 7YY
United Kingdom
Tel: +44 (0) 1638 663338
Fax: +44 (0) 1638 668890
Email: sales@nitritex.com





BioClean DB™

STERILE POUCH STYLE MASK WITH NECK GUARD

Manufacturer:	Nitritex Ltd	
Construction:	Seam:	Ultrasonically bonded
	Nosepiece:	Fully enclosed, plastic coated, soft malleable steel
	Outer facing:	Spunbond non-woven polyester
	Filter layer:	Meltblown polyester
	Inner facing:	Spunbond non-woven polyester
	Spacer:	Non-latex band
	Headbands:	Tubular knitted polyurethane Spandex yarn
	Neck guard:	Spunbond non-woven polyester
Colour:	White	
Size:	Mask:	305mm x 105mm
	Neck Guard:	210mm x 230mm
Protein Level:	Zero proteins - 100% latex free	
Sterilization:	Gamma sterilised	
Post Sterilisation Performance:	Particle filtration efficiency @ 0.3 micron	98.2%
	Bacteria filtration efficiency @ 3.0 micron	97.9%
	Differential pressure in mm H ₂ O/cm ²	3.57
Quality Systems:	Manufactured in a facility holding the following registration: ISO 9001:2000	
Sampling & Inspection:	In accordance with ISO 2859	
Storage Instructions:	Store in a cool dry place (5C-35C) away from direct sunlight and heat	
Lot Number:	This appears on the sealed bag and the outer carton to provide traceability	
Particulates	Processed in a cleanroom environment to ensure compatibility with ISO Class 4 operating conditions.	
Packaging	1 mask per inner bag; 20 inner bags per outer bag; 10 outer bags per carton (200 pieces)	
Product Codes	Re-order Code: BDBS-G	



"Comfort and Protection!"

