



VFM

BioClean Clearview™

STERILE LOOPED AND TIE-ON VISOR FACEMASKS


Why Compromise..

your sterile environment by sanitising goggles which could lead to product contamination?

BioClean Clearview™ sterile cleanroom visor facemasks will not compromise sterility and are easier to manage and control than goggles.

BioClean Clearview™ visor facemasks are 3 ply with a 3M Anti-fog visor which is optically clear and distortion free. User comfort is assured with a hydrophilic inner layer to absorb moisture and safety is enhanced by the hydrophobic outer layer that is liquid repellent.

Lint-free with very low particle release, BioClean Clearview™ visor facemasks are suitable for use in ISO Class 4 environments and are 210mm width for maximum coverage.



Your clear view to the future!



NITRITEX ASIA PTE LTD
11 Keng Cheow Street
#03-09 Riverside Piazza
Singapore 059608
Republic of Singapore
Tel: +65 6536 9123
Fax: +65 6438 5159
Email: nitritex@pacific.net.sg

NITRITEX CANADA LTD
7030 Woodbine Avenue
Suite 500, Markham
Ontario L3R 6G2
Canada
Tel: +1 (905) 946 9539
Fax: +1 (905) 946 8584
Email: info@nitritexcanada.com

NITRITEX LTD
Minton Enterprise Park
Oaks Drive, Newmarket
Suffolk, CB8 7YY
United Kingdom
Tel: +44 (0) 1638 663338
Fax: +44 (0) 1638 668890
Email: sales@nitritex.com



