



# BioClean-D™

## Disposable Lab Coat



Nitritex Ltd is pleased to present the **BioClean-D™** Lab Coat. This single-use disposable garment has been carefully researched, designed and constructed to provide the wearer with maximum comfort and protection.



The lab coat has been manufactured from a revolutionary new **Clean-Tough™** fabric which is lightweight, low linting, hard-wearing and compatible for use in ISO Class 4 cleanrooms and industrial environments.

### PPE CATEGORY III



CE 0120

The **BioClean-D™** Disposable Lab Coat is available in sterile and non-sterile versions and in a wide range of sizes to suit all users.

*"For Comfort and Protection!"*

NITRITEK ASIA PTE LTD  
11 Keng Cheow Street  
#03-09 Riverside Piazza  
Singapore 059608  
Republic of Singapore  
Tel: +65 6536 9123  
Fax: +65 6438 5159  
Email: nitritex@pacific.net.sg

NITRITEK CANADA LTD  
7030 Woodbine Avenue  
Suite 500, Markham  
Ontario L3R 6G2  
Canada  
Tel: +1 (905) 946 9539  
Fax: +1 (905) 946 8584  
Email: info@nitritexcanada.com

NITRITEK LTD  
Minton Enterprise Park  
Oaks Drive, Newmarket  
Suffolk, CB8 7YY  
United Kingdom  
Tel: +44 (0) 1638 663338  
Fax: +44 (0) 1638 668890  
Email: sales@nitritex.com





# BioClean-D™


## Disposable Lab Coat

- Manufacture:** Made and packed for Nitritex Limited
- Features:** Press stud fastening, open cuffs & rear vent
- Colour:** White
- Material:** Clean-Tough™ fabric
- Quality Systems:** Manufactured in a facility holding the following registration:  
ISO 9001:2000
- Storage Instructions:** Store in a cool dry place away from direct sunlight and heat (5°C-35°C)
- Sterilisation Method:** Where applicable - Gamma Irradiation
- Lot Number:** This appears on each individual bag and the outer carton to provide traceability
- Bar Codes:** A unique bar code (ITF 14) on the outer carton allows automatic entry into Stock Control Systems
- Packaging:** One lab coat per sealed plastic bag;  
30 plastic bags per sealed carton liner;  
One carton liner per carton. (30 lab coats)
- Printing:** Our special print encapsulation process ensures bags can be wiped down with cleaning or sterilising fluids without fear of cross contamination from printing inks.
- Product Codes:** Each size is given an individual code for re-order purposes as follows:

| SIZE | CHEST   |           | HEIGHT  |              | CODE        | CODE       |
|------|---------|-----------|---------|--------------|-------------|------------|
|      |         |           |         |              | NON-STERILE | STERILE    |
| S    | 84-92   | (33"-36") | 164-170 | (5'4"- 5'6") | BDLC S      | S-BDLC S   |
| M    | 92-100  | (36"-39") | 170-176 | (5'6"- 5'9") | BDLC M      | S-BDLC M   |
| L    | 100-108 | (39"-42") | 176-182 | (5'9"- 6'0") | BDLC L      | S-BDLC L   |
| XL   | 108-116 | (42"-45") | 182-188 | (6'0"- 6'2") | BDLC XL     | S-BDLC XL  |
| XXL  | 116-124 | (45"-48") | 188-194 | (6'2"- 6'4") | BDLC XXL    | S-BDLC XXL |



GB04/61092



*"For Comfort and Protection!"*