


YELLOW DISPOSABLE CHEMO GOWN



The yellow disposable CRG Chemo Gown is manufactured from an impervious poly-coated material and offers proven protection against cytotoxic drugs.


The ample design with waist tie ensures a comfortable gown with the barrier performance you need.

As recommended by the Workcover Guidelines for Handling Cytotoxic Drugs and Related Waste in Healthcare Establishments (NSW 1995), maximum protection is paramount for the operator's safety. The documented risks and side effects from working with cytotoxic drugs are extremely serious. The CRG chemo gown provides this protection.

Cytotoxic Impermeability: Effective barrier protection Test Results

CHEMICALS	CONCENTRATION	BREAKTHROUGH (UP TO 4 HOURS)
Etoposide	20mg/ml	None
Dacarbazine	10mg/ml	None
Vincristine Sulfate	1.0mg/ml	None

Features	Benefits
Impervious Fabric.	Suitable for use with cytotoxic and other hazardous drugs.
Elasticised cuffs.	Comfortable fit, improves ease and speed of donning gloves, hold gown in place.
Waist tie.	Easily fastens at the back for a comfortable fit.
One Size	Enables bulk buying - providing economies of scale, and eliminates sizing issues.
Low Price	Save money, less risk to change over, and excellent value for money.



Item Code: PA64001Y