

Sterile SatPax[®] 1200

Presaturated Cleanroom Wiper

COMPOSITION AND ATTRIBUTES

- Higher saturation level to support the critical solution contact time needed in the regulated industries
- 100% continuous filament polyester with knife cut edges
- Laundered and packaged in Berkshire's ISO Class 4 (Class 10) cleanroom prior to sterilization
- Gamma irradiated and validated sterile to a Sterility Assurance Level of 10^{-6} in accordance with ANSI/AAMI/ISO guidelines
- Saturated with 70% USP grade isopropyl alcohol (IPA) filtered to $0.2\mu\text{m}$
- Peel and reseal pouches allow for single wiper use and prevent evaporation
- Eliminate spray and squirt bottles and uncontrolled saturation levels
- Meeting cGMP requirements for traceability, a Certificate of Sterility is provided with each case and available on the web 24/7. Details include Date of Manufacture, Radiation Dosage, and Lot Information.

APPLICATIONS

- Compatible with ISO Class 4 (Class 10) cleanroom applications and above
- Ideal for aseptic environments and support areas where microbial exclusion is demanded
- General purpose cleaning and surface preparation
- Cleaning of exterior equipment in production areas
- Final wipe down after disinfection
- Wiping and prepping of glassware and vials

VITAL STATISTICS (AVERAGE DATA IN UNSATURATED STATE)

Basis Weight (g/m^2)	143	Extractables:	
Caliper (μm)	478	H ₂ O (%)	0.015
Fibers ($/\text{cm}^2$)	0.72	2-Propanol (%)	0.039
Particles (0.5 μm) (biaxial shake) ($\times 10^3/\text{cm}^2$)	1.6	Ions (ppm):	
Absorbency (cc/m^2)	NA	Na ⁺	NA
Specific Absorbency (cc/g)	NA	K ⁺	NA
Time to Sorption (sec)	NA	Ca ⁺⁺	NA
		Mg ⁺⁺	NA

ORDER INFORMATION

Product	Item Number	Size/Metric	Shts/Pk	Pk/Cs
Sterile SatPax 1200	SSP1200.003.12	9x9"/23x23cm	30	12